

# Eulerian Video Magnification (EVM) in Carotid Artery Image Analysis

---

Biao Jiang

March 29, 2018

Umeå University

- ① EVM Method
- ② Time-varying Texture Features
- ③ Simulations
- ④ EVM Results
- ⑤ EVM Applications

# EVM Method

---


- To magnify and perceive subtle temporal variations in videos
- Main techniques:
  - Image pyramid processing.
  - Temporal band-pass filtering to extract subtle changes.
  - Overlay the magnified image at a particular pyramid level to the original one.



# Time-varying Texture Features

---


# Texture Features



**Figure 1:** GSM colormap

- Gray-scale median (GSM)
- GLCM (Gray-level co-occurrence matrix)


# Ultrasound Intensity Cyclic Variations



**Figure 2:** Intensity cyclic variations



# GSM Variations



**Figure 3:** GSM Variations

- GSM is higher at systole than diastole.

# Simulations

---

# Effects of the Plaque Out-of-plane Motions




Figure 4

# Effects of the Plaque Compressions






Figure 5

## **EVM Results**



---

## Original Sequences vs. EVM's

# Original Frame vs. EVM's



(a)



(b)

**Figure 6:** Carotid artery B-mode images of original vs. EVM.

# **EVM Applications**

---



- Magnify the whole or parts of the plaque so that it is easier to perceive/visualize.
- Improve the segmentation performance. Since some regions in the B-mode image are magnified, i.e. image signal-to-noise (SNR) is increased, thus, for edge detection, less false contour boundaries will be produced.
- More clearly observe the differences between the healthy and stiffer arteries, between the stable and unstable plaques.
- New risk-markers. By using the texture feature analyzing methods (GSM, gray-level co-occurrence matrix (GLCM), etc.) on the post-EVM image sequences.

Thank You!